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**For simple, fast and highly accurate characterization of porous materials**

**Microtrac's BELSORP MINI X's outstanding features result in world class repeatability, with significantly reduced measurement time. The MINI X is equipped with up to 4 sample measurement ports and new high-throughput functions including multi-device control.**

The BELSORP MINI X is meticulously designed to measure specific surface area, pore size distribution and pore volume with extreme accuracy. Equipped with dedicated pressure sensors on each sample measuring port and a dedicated port for saturated vapor pressure, it enables completely independent simultaneous measurements. The new measurement software improves user productivity by streamlining measurement procedures with in software step-by-step instructions, displaying the measurement progress and sending the measurement results via e-mail. Depending on your sample throughput, a 3-port and 4-port models are offered. Additionally, the new analysis software (BELMASTER™) enables structural evaluation of a wider range of materials than ever before.

## WORLD CLASS LEVELS OF REPRODUCIBILITY

- | 多达 4 个独立测量端口和 1 个用于饱和蒸汽压测量的专用端口
- | 紧凑设计、轻量化结构
- | 专利 AFSM™ 具有无与伦比的重现性
- | Gas adsorption isotherm & NET adsorption measurement through AFSM™2 without the need of He-gas
- | Multiple Measurement Methods: High-precision measurements, multi-sample measurement, and quick BET measurement modes to suit your purpose
- | Dedicated pressure transducer for each port
- | Simultaneous control of up to 20 measurement ports via multi-device control (5 units)
- | Short measurement time with GDO
- | IoT: Process monitoring via e-mail notification system
- | 通过分子探针法进行可选的微孔分析
- | 可选的 FDA 21 CFR Part 11 合规性

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BET specific surface areas from  $0.01 \text{ m}^2/\text{g} \sim (\text{N}_2)$  and pore size distribution from  $0.7 \sim 500 \text{ nm}$  (option:  $0.35 \sim 500 \text{ nm}$  by molecular probe method).

### Simultaneous Measurement of Multiple Samples

Up to 4 specimens can be measured simultaneously, enabling the measurement times for multiple samples to be shortened significantly.

### Multiple Measurement Options


This instrument brings together highly accurate measurements with several high throughput modes (e.g. Quick BET mode, Gas Dosing Optimization). It can be used flexibly depending on customer needs, while integrated step-by-step instructions make it easy to use even for inexperienced users.

## Ease of Use: Simple Setting of Measurement Software

The BELSORP MINI X is a fully automatic measurement system that enables the user to easily set the measurement conditions using Simple Mode. This mode allows for accurate measurements with minimal sample information (e.g. pre-treatment conditions and measuring range). This is particularly useful for unknown materials. Detailed measurement configurations can be set by experienced users in Professional mode. This system allows anyone to easily obtain accurate measurement results.

## Gas Dosing Optimization (GDO)

Simple Mode includes the Gas Dosing Optimization (GDO) option, which automatically calculates the optimal gas dosing amount from previous measurement results. With this function, the measurement set points can be reliably obtained and the measurement times can be reduced significantly.

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
## Innovative Free Space Measurement for Highest Accuracy (AFSM™)

With AFSM™, the level control of liquid coolant (e. g. liquid nitrogen or argon) is not necessary anymore. Instead, we have adopted our groundbreaking, continuous free space measurement method, the patented Advanced Free Space Measurement (AFSM™). This method has achieved highest reproducibility by using a reference cell (empty sample cell; same type of sample cell used for measurement) to follow the second-to-second changes of free space. Our instrument determines an initial free space of the sample cell and reference cell. Since the change of the free space is equal in both cells (same environmental conditions), the free space change can be continuously followed by the pressure change of reference cell. Free space fluctuations caused by environmental factors can be considered:

- | LN<sub>2</sub> level changes
- | Atmospheric temperature / pressure changes
- | Temperature changes of coolant due to oxygen dissolution


With this superior technique, environmental factors can be taken into account that were not considered before.

For more information, please use our knowledge base.

### **Adsorption Measurement Without Using Helium Gas (AFSM<sup>TM2</sup>)**

In a first step the free space measurement for a blank sample cell and a blank reference cell is performed with the adsorptive gas at the appropriate adsorption temperature. Then, the adsorption measurement can be achieved by utilizing both cells without using helium gas. In the same measurement condition, since the free space changes for both cells are same, the measurement results can be obtained in shorter measurement time. There is no need for time consuming operations such as keeping the liquid level of coolant constant between the blank measurement and adsorption measurement of the free space.

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## Compact & Lightweight

Through complete optimization and rationalization of the materials, we succeeded in reducing the device size and weight.

## Optional Pretreatment of Samples at Measurement Port

For an accurate adsorption measurement, a pretreatment of the samples is required. The pretreatment process (often called an activation process) is usually performed by applying vacuum under heat to remove adsorbed gas and/or water molecules from the surface of the material without affecting the sample structure (avoid denaturation). Microtrac offers two options for sample pretreatment. Firstly, it can be performed externally using BELPREP VAC II or VAC III, which is usually preferred to increase sample throughput. Alternatively, the activation process can be performed directly at the measurement port of the BELSORP MINI X using a heater (see accessories list). Using the heater accessory, a transfer from the external pretreatment device to the measuring port can be avoided, which is an important option for sensitive (e.g. hydrophilic) samples.

## Data Analysis

The analysis software BELMASTER allows users to obtain a wide variety of analysis results through operations, such as displaying adsorption / desorption isotherms, evaluating specific surface areas using BET method etc., evaluating pore volumes (from ultra-micro- to meso- and macropores) by applying the t-plot or  $\alpha S$  method. Further, conducting mesopore analysis using the DH and BJH methods, micropore analysis using the HK or SF methods, or optional GCMC / NLDFT analysis are possible. Even more analysis options are given in our BELMASTER analysis software.

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